

INFORMATION FOR PATIENTS

Information sheet on how to use an adhesive remover spray

1. Remove packaging and protective lid from the canister. On a new canister you may have to break the seal by pressing the nozzle once. Do this facing away from the stoma
2. Lift the small tab on the stoma baseplate/flange. This tab should not be stuck to the skin. If there is no tab, gently lift the top corner of the baseplate/flange to expose the skin below
3. Hold the canister about 6-10cm away from the skin and spray the adhesive remover **between** the baseplate/flange and the skin (Figure 1)
4. Spray for 1-2 seconds. **Do not spray the adhesive remover on the top of the baseplate/flange**
5. Allow the product to drip into the space between the skin and baseplate/flange and then gently pull the baseplate/flange downwards (Figure 2)
6. When resistance is felt, spray the adhesive remover again in between the baseplate/flange and the skin for 1-2 seconds
7. Repeat this process until the baseplate/flange is completely removed
8. On average one 50ml canister of spray should last between 10-12 days if baseplate/flange changed daily.



Fig. 1



Fig. 2